



The connector has a different color depending on the diameter of the catheter.

Adhesive tab on back for reliable hanging.

Not manufactured with PVC, phthalates or latex.



LoFric catheters are coated with the unique Urotonic™ Surface Technology. This ensures the catheter remains smooth and safe, both during insertion and withdrawal.^{1,2}

Product description

LoFric® Dila-Cath™, single use, hydrophilic dilatation catheter, intended for intermittent dilatation of the urethra. The catheter does not have drainage eyes and, therefore, cannot be used for bladder emptying.

Sterilization and shelf life

Sterilized by ethylene oxide gas. Use before expiry date on package.

Package

30 catheters per customer box.
2 customer boxes per shipping box.

Material composition

Hydrophilic surface coating: Polyvinyl pyrrolidone (PVP) and sodium chloride (NaCl).

Catheter: Polyolefin-based elastomer (POBE).

Item packaging: Lacquered paper, laminated film of Polypropylene (PP)/Polyamide (PA)/Polyethylene (PE).

Product range

Reference	Type	Length	Size
406 16 XX*	LoFric® Dila-Cath™ Nelaton	40 cm	CH 16
406 18 XX*	LoFric® Dila-Cath™ Nelaton	40 cm	CH 18

*Generic article number without the two digit suffix specific for a region or country destination when distributing an article.



User-friendly packaging.



References

1. Waller et al., Spinal Cord, 1997
2. Dunger et al., Acta Endocrinologica, 1988

At Wellspect we develop innovative continence care solutions that improve quality of life for people with bladder and bowel problems. We inspire our users to build self-confidence and independence as well as good health and well-being. We have been leading the industry for over 40 years with our product brands LoFric® and Navina™. We always aim to minimize the environmental impact of our products and passionately strive to become climate neutral. We work together with users and healthcare professionals to improve clinical outcome in a sustainable way, now and for the future.

Wellspect. A Real Difference.

For more information about our products and services, please visit [Wellspect.com](https://www.wellspect.com).

Join the conversation on LinkedIn, Twitter, Facebook and Instagram.

[wellspect.com](https://www.wellspect.com)

Wellspect HealthCare, Aminogatan 1, P.O. Box 14, SE-431 21 Mölndal, Sweden.
Phone: +46 31 376 40 00. www.wellspect.com

LoFric



Manufacturer: Wellspect HealthCare, Aminogatan 1, P.O. Box 14, SE-431 21 Mölndal, Sweden.
Phone: +46 31 376 40 00. www.wellspect.com

CE 2797

Rx Only