

Navina Classic Product Information

Product description

Transanal Irrigation System with manual control unit.

Shelf life and uses

All Navina Irrigation System components have a shelf life of 3.5 years, except the Navina Catheters and the Navina Cone, which have 2 years. When Navina components are packed in sets, the component with shortest shelf life determine the shelf life of the set.

Navina catheters and cones are single use, the water container has 15 uses, and the Navina Classic control unit and tube set has 100 uses.

Material composition

Catheter Material: Styrene Ethylene Butylene Styrene (SEBS)/ Polypropylene (PP) with a hydrophilic Polyvinylpyrrolidone (PVP) coating.

Primary Package Material: PE/PET Pouch.

Secondary Package Material: Cardboard and corrugated cardboard.

Intended use

The Navina Systems is intended for Transanal Irrigation by instilling water up into the lower part of the colon through a rectal catheter.

Indications for use

The Navina Systems is indicated to help adults and children from 3 years who suffer from fecal incontinence, chronic constipation and/or time consuming bowel management. By instilling water up into the lower part of the colon, the peristaltic muscles in the bowel can be triggered and start to evacuate the content of the lower colon and rectum.

Contraindications

Do NOT use Navina Systems if you have one or more of the following:

- Known anal or colorectal stenosis.
- Active inflammatory bowel disease.
- Acute diverticulitis.
- Colorectal cancer.
- Ischemic colitis.
- You are within three months of anal or colorectal surgery.
- You are within 4 weeks of previous endoscopic polypectomy.
- You are pregnant.

As the list may not be exhaustive, healthcare professionals will always consider individual user factors as well.



Navina Classic System

Reference	Type	Specification
69005 69007 69017	Navina Classic System Regular Navina Classic System Small Navina Classic System Cone	1 Navina Classic control unit 1 Water container 1 Tube set 2 Rectal catheters 1 Navina case 1 Accessory set
69006 69008 69018	Navina Classic System Regular Refill Navina Classic System Small Refill Navina Classic System Cone Refill	1 Navina Classic control unit 1 Water container 1 Tube set 2 Rectal catheters
69003 69004 69016	Navina Consumable Set Regular Navina Consumable Set Small Navina Consumable Set Cone	1 Water container 15 Rectal catheters
68940 68943	Navina Catheter Set Regular Navina Catheter Set Small	10 Rectal catheters
69011	Navina Tube Set	1 Catheter tube 1 Water container tube
69012	Navina Accessories	2 Grip rings (1 vertical, 1 horizontal) 1 Lanyard 1 Positioning strap 1 Navina Smart positioning clip
69013	Navina Classic control unit	1 Navina Classic control unit

Navina™ Classic



Manufacturer: Wellspect HealthCare, Aminogatan 1, P.O. Box 14, SE-431 21 Mölndal, Sweden.
Phone: +46 31 376 40 00. www.wellspect.com



Rx Only