



*Discreet design.*

*User-friendly handle for good grip and non-touch technique.*

*Adhesive tab on back for reliable hanging.*

*Instant activation - sterile salt solution included, squeeze and it's ready to use.*

*Peel-open tab.*

*The outer packaging doubles as a hygienic and discreet disposal pouch after use.*

*Not manufactured with PVC, phthalates or latex.*

*Smooth eyelets.<sup>1,2</sup>*



*LoFric catheters are coated with the unique Urotonic™ Surface Technology. This ensures the catheter remains smooth and safe, both during insertion and withdrawal.<sup>3,4</sup>*

### Product description

LoFric® Sense™, single use, hydrophilic urinary catheter with salt solution for activation, intended for intermittent urinary catheterization.

### Sterilization and shelf life

Sterilized by irradiation. Use before expiry date on package.

### Package

30 catheters per customer box. 4 customer boxes per shipping box.  
5 catheters per sample pack.

### Material composition

**Hydrophilic surface coating:** Polyvinyl pyrrolidone (PVP).

**Catheter:** Polyolefin-based elastomer (POBE).

**Sachet:** Laminate foil of polyethylene terephthalate (PET), aluminum and polyethylene (PE). Sterile water with sodium chloride (NaCl).

**Item packaging:** Laminate foil of polyethylene terephthalate (PET) and polyethylene (PE).

### Product range

Reference	Type	Length	Size
416 08 XX*	LoFric® Sense™ Nelaton (female)	15 cm	CH 08
416 10 XX*	LoFric® Sense™ Nelaton (female)	15 cm	CH 10
416 12 XX*	LoFric® Sense™ Nelaton (female)	15 cm	CH 12
416 14 XX*	LoFric® Sense™ Nelaton (female)	15 cm	CH 14
416 08 05	LoFric® Sense™ Nelaton (female)	15 cm	CH 08 <small>sample pack</small>
416 10 05	LoFric® Sense™ Nelaton (female)	15 cm	CH 10 <small>sample pack</small>
416 12 05	LoFric® Sense™ Nelaton (female)	15 cm	CH 12 <small>sample pack</small>
416 14 05	LoFric® Sense™ Nelaton (female)	15 cm	CH 14 <small>sample pack</small>

*\*Generic article number without the two digit suffix specific for a region or country destination when distributing an article.*

*The effective length is 10 cm. Effective length, length of the catheter that can be inserted into the body.*



*User-friendly packaging.*



*Sample pack, 5 catheters.*



**LoFric® Sense™**

## References

1. Lundgren et al., Spinal Cord, 2000
2. Newman et al., Springer, 2018
3. Waller et al., Spinal Cord, 1997
4. Dunger et al., Acta Endocrinologica, 1988

At Wellspect we develop innovative continence care solutions that improve quality of life for people with bladder and bowel problems. We inspire our users to build self-confidence and independence as well as good health and well-being. We have been leading the industry for over 40 years with our product brands LoFric® and Navina™. We always aim to minimize the environmental impact of our products and passionately strive to become climate neutral. We work together with users and healthcare professionals to improve clinical outcome in a sustainable way, now and for the future.

### **Wellspect. A Real Difference.**

For more information about our products and services, please visit [Wellspect.com](https://www.wellspect.com).

Join the conversation on LinkedIn, Twitter, Facebook and Instagram.

**wellspect.com**

Wellspect HealthCare, Aminogatan 1, P.O. Box 14, SE-431 21 Mölndal, Sweden.  
Phone: +46 31 376 40 00. [www.wellspect.com](https://www.wellspect.com)

LoFric



Manufacturer: Wellspect HealthCare, Aminogatan 1, P.O. Box 14, SE-431 21 Mölndal, Sweden.  
Phone: +46 31 376 40 00. [www.wellspect.com](https://www.wellspect.com)

CE 2797

Rx Only