



Luer-lock. A connector type that is safe and minimizes the risk of spillage and unintentional exposure. Instillation solution is held completely enclosed.**

The connector has a different color depending on the diameter of the catheter.

Adhesive tab on back for reliable hanging.

Not manufactured with PVC, phthalates or latex.



LoFric catheters are coated with the unique Urotonic™ Surface Technology. This ensures the catheter remains smooth and safe, both during insertion and withdrawal.^{1,2}

Smooth eyelets.^{3,4}



Nelaton Tiemann

Product description

LoFric® Insti-Cath™, single use, hydrophilic urinary catheter, intended for intermittent bladder instillation therapy.

Sterilization and shelf life

Sterilized by ethylene oxide gas. Use before expiry date on package.

Package

30 or 10 catheters per customer box. 4 or 7 customer boxes per shipping box.

Material composition

Hydrophilic surface coating: Polyvinyl pyrrolidone (PVP) and sodium chloride (NaCl).

Catheter: Polyolefin-based elastomer (POBE).

Item packaging: Lacquered paper, laminated film of Polypropylene (PP) / Polyamide (PA) / Polyethylene (PE).

Product range

Reference	Type	Length	Size
407 08 XX*	LoFric® Insti-Cath™ Nelaton (male)**	40 cm	CH 08
407 10 XX*	LoFric® Insti-Cath™ Nelaton (male)**	40 cm	CH 10
407 12 XX*	LoFric® Insti-Cath™ Nelaton (male)**	40 cm	CH 12
407 14 XX*	LoFric® Insti-Cath™ Nelaton (male)**	40 cm	CH 14
408 08 XX*	LoFric® Insti-Cath™ Nelaton (female)**	20 cm	CH 08
408 10 XX*	LoFric® Insti-Cath™ Nelaton (female)**	20 cm	CH 10
408 12 XX*	LoFric® Insti-Cath™ Nelaton (female)**	20 cm	CH 12
408 14 XX*	LoFric® Insti-Cath™ Nelaton (female)**	20 cm	CH 14
409 12 XX*	LoFric® Insti-Cath™ Tiemann (male)**	40 cm	CH 12
440 08 XX*	LoFric® Insti-Cath™ Nelaton (paediatric)**	20 cm	CH 08
440 10 XX*	LoFric® Insti-Cath™ Nelaton (paediatric)**	20 cm	CH 10

*Generic article number without the two digit suffix specific for a region or country destination when distributing an article.

**Only instill solutions approved for intravesical bladder administration and according to instructions.



User-friendly packaging.

LoFric® Insti-Cath™

References

1. Waller et al., Spinal Cord, 1997
2. Dunger et al., Acta Endocrinologica, 1988
3. Lundgren et al., Spinal Cord, 2000
4. Newman et al., Springer, 2018

At Wellspect we develop innovative continence care solutions that improve quality of life for people with bladder and bowel problems. We inspire our users to build self-confidence and independence as well as good health and well-being. We have been leading the industry for over 40 years with our product brands LoFric® and Navina™. We always aim to minimize the environmental impact of our products and passionately strive to become climate neutral. We work together with users and healthcare professionals to improve clinical outcome in a sustainable way, now and for the future.

Wellspect. A Real Difference.

For more information about our products and services, please visit [Wellspect.com](https://www.wellspect.com).

Join the conversation on LinkedIn, Twitter, Facebook and Instagram.

[wellspect.com](https://www.wellspect.com)

Wellspect HealthCare, Aminogatan 1, P.O. Box 14, SE-431 21 Mölndal, Sweden.
Phone: +46 31 376 40 00. www.wellspect.com

LoFric



Manufacturer: Wellspect HealthCare, Aminogatan 1, P.O. Box 14, SE-431 21 Mölndal, Sweden.
Phone: +46 31 376 40 00. www.wellspect.com

CE 2797

Rx Only